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(一) US EPA, Colorimetric Screening Method for RDX and HMX in

Soil, SW-846 Method 8510, 2007.

- (二) M. E. Walsh and T. F. Jenkins, Development of a Field Screening Method for RDX in Soil, US Army Cold Regions Research and Engineering Laboratory, Special Report 91-7, Hanover, New Hampshire 03755, 1991.
- (三) T. F. Jenkins, M. E. Walsh, P. W. Schumacher and P. G. Thorne, Development of Colorimetric Field Screening Methods for Munitions Compounds in Soil, Society of Photo-Optical Instrumentation Engineers, 2504, 324 - 333, 1995.
- (四) RDX Soil Test System Users Guide, EnSys, Inc., 1996.
- (五) S. P. Arrowood, P. P. McDonald and J. P. Mapes, RDX Soil Test Kit Validation Results, EnSys, Inc., March, 1993.
- (六) Onsite Analytical Technologies Umatilla Army Depot," USEPA Project No. 71370, Prepared by Black and Veatch Special Projects Corporation for USEPA Region 10, Draft May, 1997.
- (七) T. F. Jenkins, M. E. Walsh, P. G. Thorne, S. Thiboutot, G. Ampleman, T. A. Ranney and C. L. Grant, Assessment of Sampling Error Associated with Collection and Analysis of Soil Samples at a Firing Range Contaminated with HMX, US Army Cold Regions Research and Engineering Laboratory, Special Report 97-22, Hanover, New Hampshire 03755, September, 1997.

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